

## Product overview


# Sialidase Testing Panel – Procainamide Labelled Glycan Standard

Cat. # CPROC-STP-NEUAC-01

### Description


Sialidase Testing Panel is a glycan standard containing a mixture of  $\alpha$ 2-3,  $\alpha$ 2-6 and  $\alpha$ 2-8 sialylated oligosaccharides:

1. 3'-Sialyl Lewis X [**Neu5Ac-a2-3Gal-b1-4(Fuc-a1-3)GlcNAc**] – contains  $\alpha$ 2-3 linked sialic acid along with branched  $\alpha$ 1-3 fucose (introducing steric hindrance for some enzymatic reactions)
2. Disialyllactose [**Neu5Ac-a2-8NeuAc-a2-3Gal-b1-4Glc**] – linear oligosaccharide containing  $\alpha$ 2-8 linked sialic acid
3. Sialyllacto-N-tetraose a [**Neu5Ac-a2-3Gal- $\beta$ 1-3GlcNAc- $\beta$ 1-3Gal- $\beta$ 1-4Glc**] – linear oligosaccharide containing  $\alpha$ 2-3 sialic acid attached to galactose
4. Sialyllacto-N-tetraose c [**Neu5Ac-a2-6Gal-b1-4GlcNAc-b1-3Gal-b1-4Glc**] – linear oligosaccharide containing  $\alpha$ 2-6 sialic acid attached to galactose
5. Disialyllacto-N-tetraose [**Neu5Ac-a2-3Gal-b1-3(Neu5Ac-a2-6)GlcNAc-b1-3Gal-b1-4Glc**] – branched sialic acid oligosaccharide




**Figure 1:** HILIC-UPLC profiles of procainamide labelled Sialidase Testing Panel (CPROC-STP-NEUAC-01) in comparison with GHP (CPROC-GHP-30).

|                         |   |
|-------------------------|---|
| <b>Form</b>             | Dry. Dried by centrifugal evaporation from an aqueous solution.   |
| <b>Size</b>             | 20 pmol   |
| <b>Molecular weight</b> | <ol style="list-style-type: none"> <li>1. 3'-Sialyl Lewis X 1039.4661</li> <li>2. Disialyllactose 1143.4770</li> <li>3. Sialyllacto-N-tetraose a 1217.5138</li> <li>4. Sialyllacto-N-tetraose c 1217.5138</li> <li>5. Disialyllacto-N-tetraose 1508.6092</li> </ol>   |
| <b>Purity</b>           | 89% pure as assessed by UHPLC.  |
| <b>Storage</b>          | Store at -20°C both before and after dissolution. This product is stable for at least 5 years as supplied.  |
| <b>Shipping</b>         | The product can be shipped at ambient when dry. After dissolution, ship on dry ice.   |
| <b>Handling</b>         | Allow the unopened vial to reach ambient temperature. Centrifuge to ensure that most of the lyophilized material is at the bottom of the vial. Add the desired volume of reconstitution medium, re-cap and mix thoroughly to bring all the oligosaccharide into solution. Ensure that any glass, plasticware or solvents used are free of glycosidases and environmental carbohydrates. |
| <b>Example use</b>      | Sialidase Testing Panel can be used as process positive control for sialidase digestions. Applied alongside the samples, it enables to test if the sialidase used has required specificity and if it had worked correctly.  |



**Figure 2:** Workflow for exoglycosidase sequencing of glycans based on HPLC/LC-MS detection. Sialidase Testing Panel can be incorporated as a process control for sialidase digestion. Other controls e.g. fucosylated,  $\alpha$ -galactosylated,  $\beta$ -galactosylated glycan standards can be incorporated alongside, depending on exoglycosidases used. Example methods of detection for reaction products include MS and FLD-LC.



**Figure 3:** HILIC-UPLC profiles of procainamide labelled Sialidase Testing Panel (CPROC-STP-NEUAC-01). Top chromatogram shows undigested sample, bottom chromatogram shows sample treated with broad specificity sialidase (E-S001). Peaks were labelled with glycan structures. Arrows illustrate glycan digestion pathways.

## Related products

| Ludger Cat. No.  | Description                            |
|------------------|--|
| CAB-STP-NEUAC-01 | Sialidase Testing Panel, 2-AB labelled |
| E-S001           | Sialidase Au $\alpha$ (2-3,6,8,9)      |
| E-S005           | Sialidase Cp $\alpha$ (2-3,6)          |
| E-S007           | Sialidase Sp $\alpha$ (2-3)            |

## Warranties and liabilities

Ludger warrants that the above product conforms to the attached analytical documents. Should the product fail for reasons other than through misuse Ludger will, at its option, replace free of charge or refund the purchase price. This warranty is exclusive and Ludger makes no other warranties, expressed or implied, including any implied conditions or warranties of merchantability or fitness for any particular purpose.

Ludger shall not be liable for any incidental, consequential or contingent damages.

This product is intended for *in vitro* research only. Not for human or drug use.

## Document Revision Number

Document # CPROC-STP-NEUAC-01-Product-Overview-v1.0