

Product Guide for LudgerTag™ Procainamide Glycan Labeling Kit with sodium cyanoborohydride

# **Product # LT-KPROC-96**

Procainamide Cat # LT-PROC-96 Batch # B76L-07 Store at Boom Tem

Ludger Document # LT-KPROC-96-Guide-v1.0

# Ludger Ltd

Culham Science Centre Oxford OX14 3EB United Kingdom

Tel: +44 1865 408 554 Fax: +44 870 163 4620 Email: info@ludger.com

# www.ludger.com



# Contents

# Page

Contents	2
Specifications for LT-KPROC-96	3
Kit Contents	4
Additional Reagents and Equipment Required	4
Time Line for Labeling	4
Printing Out Method Instructions For Your Lab	4
Procainamide Glycan Labeling System	5
Labeling Method	5
LudgerClean™ Post-Labeling Sample Clean-up	8
Analysis of LudgerTag™ Procainamide Labeled Glycans	9
References1	1
The Reductive Amination Reaction1	2
Warranties and Liabilities1	3
Document Revision Number1	3
SAFETY DATA SHEET	4
Version: 1.0	4
SAFETY DATA SHEET	2
Version: 2.0	2
SAFETY DATA SHEET	9
Version: 3.0	9



# **Specifications for LT-KPROC-96**

Application	For labeling of free glycans with procainamide (PROC).
Description	The kit contains reagents for the conjugation of dye to the free reducing end of the glycan by a reductive amination reaction.
Dye Properties	Mass free dye = 235.33.

Fluorescence,  $\lambda_{ex}$  = 310 nm,  $\lambda_{em}$  = 370 nm.



96 separate analytical samples per set of labeling reagents
From 25 pmol up to 25 nmol glycans per sample.
Any purified glycans with free reducing termini can be labeled.
No detectable (< 2 mole per cent) loss of sialic acid, fucose, sulfate, or phosphate.
Essentially stoichiometric labeling.
Store at room temperature in the dark. Protect from sources of heat, light, and moisture. The reagents are stable for at least two years as supplied.
The product can be shipped at ambient temperature.
Ensure that any glass, plasticware or solvents used are free of glycosidases and environmental carbohydrates. Use powder-free gloves for all sample handling procedures and avoid contamination with environmental carbohydrate.
All steps involving labeling reagents must be performed in a dry environment with dry glassware and plasticware. Once individual vials of reagents are opened, their contents should be used immediately and excess then discarded according to local safety rules.
For research use only. Not for human or drug use
Please read the Safety Data Sheets (SDS's) for all chemicals used. All processes involving labeling reagents should be performed using appropriate personal safety protection - eyeglasses, chemically resistant gloves (e.g. nitrile), etc and where appropriate in a laboratory fume cupboard.



# **Kit Contents**



Each kit contains two labeling reaction sets. Each labeling reaction set consists of one vial of each of the following:

Cat. #	ltem	Quantity
LT-PROC-96	procainamide dye	130 mg
LT-CYANOB-96	sodium cyanoborohydride	72 mg
LT-ACETIC-DMSO-96	30% acetic acid in DMSO	2000 μL

# **Additional Reagents and Equipment Required**

- Pure water •
- Heating block, oven, or similar dry heater (a water bath cannot be used) set at 65°C .

- Centrifugal evaporator (e.g. Savant, Heto, or similar) •
- Reaction vials (e.g. polypropylene microcentrifuge vials) •
- Note: Further consumables are required if doing the optional post-labeling sample cleanup (see • Section on Sample Cleanup)

# **Time Line for Labeling**

The LudgerTag<sup>™</sup> labeling procedure takes 2.5 hours with just 1 hour for the actual labeling incubation.

Procedure	Time	Elapsed Time (hours)
Transfer samples to reaction tube/PCR plate and dry	60 min	1
Make up and add labeling reagent	30 min	1.5
Incubate samples with reagent	1 hour	2.5

# **Printing Out Method Instructions For Your Lab**

If you would like to print out this method for your lab then please print out pages 5, 6 and 7. These pages contain all the information required to use the kit.



# **Procainamide Glycan Labeling System**



For more information on our Procainamide system and to view a presentation, please visit <u>www.ludger.com/procainamide</u>

# **Labeling Method**

# **1** Purify the glycans

Before labeling the glycans, it is preferable to remove non-carbohydrate contaminants from the samples such as protein/peptide material, salts and detergents. This can be achieved using LudgerClean EB10 cartridges (<u>LC-EB10-A6</u>) or a LudgerClean Protein Binding Membrane (<u>LC-PBM-96</u>). Clean up of O-glycans can be achieved using LudgerClean CEX cartridges (<u>LC-CEX-A6</u>).

# 2 Transfer samples to reaction vials or a PCR plate



The kit is designed to label up to 25 nmols of glycans per reaction. With a single pure glycan, as
little as 5 pmoles per reaction can be labeled and detected in subsequent HPLC and MS analysis.
Suitable reaction vials include small polypropylene microcentrifuge vials, vials for PCR work or 96
well PCR plates if performing higher throughput sample analysis.



# 3 Dry the samples



Drying down of the samples should be carried out using a centrifugal evaporator. If this is not possible then freeze drying (lyophilization) can be used with caution (in particular, ensure that the sample dries to a small, compact mass at the very bottom of the vial). Do not subject samples to high temperatures (>28°C) or extremes of pH as these conditions could result in acid catalysed loss of sialic acids (high temperatures, low pH) or epimerization of the glycan reducing terminus (at high pH).

# 4 Prepare the dye solution



- Add 1.2 mL of kit component LT-ACETIC-DMSO-96 (30% acetic acid in DMSO) to the bottle of procainamide dye (LT-PROC-96), replace the lid and shake and vortex until the dye is dissolved.
- Add 1.2 mL of water to the dissolved dye solution, replace the lid and shake and vortex.

• Transfer the 2.4 mL of dissolved dye solution to a vial of reductant (LT-CYANOB-96) replace the lid and shake and vortex until the reductant is dissolved.

# 🐗 Ludger

# 5 Add labeling reagent to samples



 Add 20 µL of labeling reagent to each glycan sample, cap the vials or PCR plate, mix thoroughly and centrifuge briefly (5-10 seconds) to ensure the labeling solution is at the bottom of the vial.





# 6 Incubate



 Place the reaction vials in an oven or heating block set at 65°C and incubate for 1 hour. We recommend using an oven for the incubation step. An oven provides all around heat for the reaction vials.

The samples must be completely dissolved in the labeling solution for efficient labeling. To encourage complete dissolution the samples can be vortexed 15 minutes after the start of the incubation then continued.

# 7 Cool and Centrifuge



• After the incubation period remove the samples from the oven or heating block, allow them to cool to room temperature then briefly centrifuge the reaction vials to ensure the sample is at the bottom of the vial and not in the lid.



# LudgerClean<sup>™</sup> Post-Labeling Sample Clean-up



For most applications, we recommend that you perform post-labeling sample cleanup to remove nonglycan material, such as excess dye and other labeling reagents, prior to analysis by HPLC, MS or LC-MS.

# Benefits of post-labeling cleanup

- Removes free dye and chemicals that may interfere with sample analysis by HPLC and LC-MS.
- Lengthens lifespan of HPLC column.
- Smaller glycans such as O-links close to the start of an HPLC gradient will be detected without interference from the free dye.
- A full range of HPLC columns can be used.
- Sample capacity issues of HPLC columns less likely.

Post-labeling sample clean-up of procainamide (PROC) labeled N-glycans can be achieved using either LudgerClean<sup>™</sup> S cartridges (LC-S-A6) or a LudgerClean Procainamide Clean-up Plate (LC-PROC-96) if performing higher throughput sample analysis.

Post-labeling sample clean-up of PROC labeled O-glycans can be achieved using LudgerClean<sup>™</sup> S cartridges (LC-S-A6).

This is summarised in Table 1.

Type of glycan	LudgerClean Product	Product Code
PROC labeled N-Glycans	LudgerClean Procainamide Clean-up Plate	LC-PROC-96
	LudgerClean™ S cartridges	<u>LC-S-A6</u>
PROC labeled O-Glycans	LudgerClean™ S cartridges	LC-S-A6

Table 1: Selection of LudgerClean products for clean-up of PROC labeled N- or O-glycans

# Analysis of LudgerTag<sup>™</sup> Procainamide Labeled Glycans

LudgerTag<sup>™</sup> procainamide labeled glycans may be studied by a number of different analytical methods including (U)HPLC, ESI mass spectrometry and LC-ESI-MS.

Example UHPLC conditions – conditions vary from instrument to instrument. Check these conditions are optimal for your equipment.

- Sample made up in solvent equivalent to starting gradient e.g. 76% acetonitrile.
- Column: Waters BEH Glycan column 15 cm x 2.1 mm.
- Column temperature: 40°C
- Fluorescence detector settings: Excitation wavelength: 310 nm, Emission wavelength: 370 nm.
- Solvent A: 50 mM ammonium formate buffer pH 4.4 (Ludger Product: LS-N-BUFFX40)
- Solvent B: Acetonitrile

Time (min)	% solvent B	Flow Rate (mL/min)
0.00	76	0.4
53.5	51	0.4
55.5	0	0.25
57.5	0	0.25
59.5	76	0.25
65.5	76	0.25
66.5	76	0.4
70.0	76	0.4

 Table 2: Long UHPLC gradient for samples where the glycan profile is unknown

 and where glycans may be large and/or highly sialylated.

The use of this UHPLC gradient is demonstrated for the separation of N-glycans that have been released from human IgG and erythropoietin followed by procainamide labeling (LT-KPROC-96) and clean-up (LC-PROC-96) (Figures 1 and 2). Once the range of glycans present in your sample has been determined the UHPLC gradient can be optimised and significantly shortened (down to 15 minutes on some UHPLC instruments).

LC-MS analysis of procainamide labeled N-glycans allows the user to obtain both a fluorescent chromatogram of the separated N-glycans (Figure 1) and an ESI-MS spectrum that can be combined to give a greater overall picture of N-glycan structures present. The MS data can be used to assign m/z masses to each of the peaks in the fluorescent chromatogram thereby giving possible structural identity, in the form of monosaccharide composition, to each of the peaks. Single peaks in the fluorescent chromatogram can be determined to have multiple glycan structures present by ESI-MS analysis.





Figure 1: Procainamide labeled human IgG glycans after LC-PROC-96 clean-up



Figure 2: Procainamide labeled erythropoietin glycans after LC-PROC-96 clean-up



# References

1. Kozak RP, Tortosa CB, Fernandes DL, Spencer DI. Comparison of procainamide and 2aminobenzamide labeling for profiling and identification of glycans by liquid chromatography with fluorescence detection coupled to electrospray ionization-mass spectrometry. *Anal Biochem.* 2015; 486: 38-40. doi: 10.1016/j.ab.2015.06.006.

2. Klapoetke S, Zhang J, Becht S, Gu X, Ding X. The evaluation of a novel approach for the profiling and identification of N-linked glycan with a procainamide tag by HPLC with fluorescent and mass spectrometric detection. *J Pharm Biomed Anal.* 2010; 53(3): 315-24. doi: 10.1016/j.jpba.2010.03.045.

3. Pabst M, Kolarich D, Pöltl G, Dalik T, Lubec G, Hofinger A, Altmann F. Comparison of fluorescent labels for oligosaccharides and introduction of a new postlabeling purification method. *Anal Biochem.* 2009; 384(2): 263-73. doi: 10.1016/j.ab.2008.09.041.

4. Liu R, Giddens J, McClung CM, Magnelli PE, Wang LX, Guthrie EP. Evaluation of a glycoengineered monoclonal antibody via LC-MS analysis in combination with multiple enzymatic digestion. *MAbs.* 2016; 8(2): 340-6. doi: 10.1080/19420862.2015.1113361.



# **The Reductive Amination Reaction**

# 1. Schiff's base formation.

This requires a glycan with a free reducing terminus which is equilibrium between the ring closed (cyclic) and ring open (acyclic) forms. The primary amino group of the dye performs a nucleophilic attack on the carbonyl carbon of the acyclic reducing terminal residue to form a partially stable Schiff's base.

# 2. Reduction of the Schiff's base.

The Schiff's base imine group is chemically reduced to give a stable labeled glycan.



Figure 3: Labeling of a glycan with procainamide (PROC) by reductive amination.



# Warranties and Liabilities

Ludger warrants that the above product conforms to the attached analytical documents. Should the product fail for reasons other than through misuse Ludger will, at its option, replace free of charge or refund the purchase price. This warranty is exclusive and Ludger makes no other warrants, expressed or implied, including any implied conditions or warranties of merchantability or fitness for any particular purpose.

Ludger shall not be liable for any incidental, consequential or contingent damages.

This product is intended for in vitro research only.

# **Document Revision Number**

Document # LT-KPROC-96-Guide-v1.0

# SAFETY DATA SHEET

Version: 1.0 Date written: 31<sup>st</sup> March 2015 Date reviewed: 6<sup>th</sup> Mar 2017

# SECTION 1. IDENTIFICATION OF THE SUBSTANCE/PREPARATION AND OF THE COMPANY / UNDERTAKING

LT-ACETIC-DMSO-96

Product Name

Acetic Acid / dimethyl sulfoxide solution

Product Catalogue Name

Company:

Telephone: Emergency Telephone: Email: Ludger Ltd Culham Science Centre Abingdon Oxford OX14 3EB 01865 408554 01865 408554 info@ludger.com

# **SECTION 2. HAZARDS IDENTIFICATION**

#### 2.1 Classification of the substance or mixture Classification according to Regulation (EC) No. 1272/2008 [EU-GHS-CLP] Flammable liquids (Category 3) Skin corrosion (Category 1A)

#### 2.2 Label elements



Signal Word: Danger

Hazard Statement(s)	
H226	Flammable liquid and vapour
H314	Causes severe skin burns and eye damage.
Precautionary Statem	ent(s)
P280 protection.	Wear proactive gloves/ protective clothing/ eye protection/ face
P305+P351+P338 contact	IF IN EYES: Rinse cautiously with water for several minutes. Remove
	lenses, if present and safe to do so. Continue rinsing.
P310	Immediately call a POISON CENTRE or doctor/ physician.

# 2.3 Other hazard information:

None

# SECTION 3. COMPOSITION/INFORMATION ON INGREDIENTS

#### 3.1 Substances

Synonyms:

DMSO, methyl sulfoxide, dimethyl sulfoxide Glacial acetic acid



Formula:	DMSC
Molecular Weight:	Acetic DMSC
-	Acetic

 $C_2H_6OS$ Acid:  $C_2H_4O_2$ D: 78.13 g/mol Acid: 60.05 g/mol

Component		Concentration	
Name	Dimethyl Sulfoxide	70%	
CAS-No.	67-68-5		
EC-No.	200-664-3		
Name	Acetic Acid	30%	
CAS-No.	64-19-7		
EC-No.	200-580-7		
Index-No.	607-002-00-6		

# **SECTION 4. FIRST AID MEASURES**

#### 4.1 Description of first aid measures

#### **General Advice**

Consult a physician if exposure causes ill effects and if in any doubt. Show this safety data sheet to the physician/ first responder in attendance.

#### If Ingested

Do NOT induce vomiting. Rinse mouth well with water. Never give anything by mouth to an unconscious person.

#### If skin is exposed

Remove all contaminated clothing immediately; wash the area well with plenty of soap and water.

#### If eyes are exposed

Flush eyes with plenty of water/ eye wash solution for at least 15 minutes, if present and safe to do so, remove contact lenses and continue rinsing.

#### If inhaled

Move effect person to fresh air. If not breathing give artificial respiration.

#### 4.2 Most important symptoms and effects, both acute and delayed

Nausea, Fatigue and Headache. To the best of our knowledge, the chemical, physical and toxicological properties have not been thoroughly investigated.

#### 4.3 Indication of immediate medical attention and special treatment needed No data available.

# **SECTION 5. FIRE-FIGHTING MEASURES**

#### 5.1 Extinguishing media

Small fires: Use extinguishing media such as "alcohol" foam, dry chemical or carbon dioxide. Large fires: Use extinguishing media such as water, from a far away distance as possible. Use very large quantities of water as mist or spray to flood the fire and the combustible material. Cool all affected containers with large quantities of water.

#### 5.2 Special hazards arising from the substance or mixture

Carbon oxides, Sulphur oxides

#### 5.3 Advice for fire fighters



Wear self contained breathing apparatus fir fire fighting if necessary, to spray cool water on any unopened containers near the fire.

# **SECTION 6. ACCIDENTAL RELEASE MEASURES**

#### 6.1 Personal precautions, protective equipment and emergency procedures

Avoid breathing vapours, gas or mist. Remove all sources of ignition. Beware of vapours accumulating to form explosive concentrations. Vapours can accumulate in low areas.

#### **6.2 Environmental Precautions**

Prevent further leakage or spillage if safe to do so, e.g. with spill mats. Do not let the product enter drains.

#### 6.3 Methods and material for containment and cleaning up

Contain the spillage and put the collected material into a suitable container with a secure lid. Wash the area well, do not let run off into the drains, collect as waste.

#### 6.4 Reference to other sections

See section 13 for disposal of waste material(s).

# **SECTION 7. HANDLING AND STORAGE**

#### 7.1 Precautions for safe handling

Avoid inhalation of vapour or mist. Keep away from sources of ignition- No smoking. Take measures to prevent the build up of electrostatic charge.

#### 7.2 Conditions for safe storage, including any incompatibilities

Store in a cool place. Keep container closed in a dry well ventilated place.

#### 7.3 Specific end uses

No data available

## **SECTION 8. EXPOSURE CONTROLS/PERSONAL PROTECTION**

#### 8.1 Control parameters

Components with workplace control parameters.

#### ACETIC ACID

CAS-No.	Value	Control Parameters	Update	Basis
64-19-7	TWA	10ppm 25mg/m3	1991-07-05	Europe. Commission Directive 91/322/EEC on establishing indicative limit on values.
Remarks	Indicative			

DMSO contains no substances with occupational exposure limit values.

#### 8.2 Exposure controls

#### Appropriate engineering controls

Handle in accordance with good laboratory hygiene and safety practice. Wash hands before breaks and at the end of the day.

#### Personal Protective Equipment Eye / face protection

Use equipment for eye protection tested and approved under appropriate government standards such as NIOSH (US) or EN 166 (EU).



#### **Skin protection**

Handle with gloves, which should be inspected before use. Use proper glove removal technique (removal without the outside of the glove touching the skin) to avoid contact with the skin/chemical. Dispose of contaminated gloves as Laboratory waste in accordance with applicable laws and good laboratory practices. Wash and dry hands.

Gloves should be of the standard that will stratify the specifications of EU directive 89/696/EEC and the standard EN 374 derived from it.

#### **Body Protection**

The type of protective clothing must be selected according to the amount of substance at the specific workplace being used. Impervious coats or laboratory coats.

#### **Respiratory protection**

Use substance in an operation fume hood/ outside venting extraction cupboard. Wear full face respirator if appropriate to use, must be tested and approved under appropriate government standards such as NIOSH (US) or CEN (EU).

## **SECTION 9. PHYSICAL AND CHEMICAL PROPERTIES**

#### 9.1 Information on basic physical and chemical properties

Appearance

Odour Odour threshold bН Freezing/Melting Point Initial boiling point and boiling range Flash Point Evaporation rate Flammability Upper/lower flammability or explosive limits Vapour Pressure, Pa at temperature degree C **Relative Density** Solubility in water and solvents Partition coefficient: n-octanol/water Auto ignition temperature Decomposition temperature Viscosity Explosive properties Oxidising properties

Form: Liquid, clear Colour: Colourless Strong No data available Completely miscible No data available No data available

# 9.2 Other information

No data available

# SECTION 10. STABILITY AND REACTIVITY

10.1 ReactivityNo data available10.2 Chemical stabilityNo data available

**10.3 Possibility of hazardous reactions** No data available

#### **10.4 Conditions to Avoid**



#### Heat, flames and sparks

#### **10.5 Incompatible materials**

Acid chlorides, Phosphorus halides, Strong oxidizing agents and strong reducing agents, soluble carbonates and phosphates, hydroxides, metals, peroxides, permanganates, e.g. potassium permanganate, Amines and Alcohols.

#### **10.6 Hazardous decomposition products**

Other decomposition products - No data available

# SECTION 11. TOXICOLOGICAL INFORMATION

11.1 Information on toxicological effects DMSO Acute toxicity

LD50 Oral – Řat – 14,500mg/kg LC50 Inhalation – Rat – 4h – 40250ppm LD50 Dermal – Rabbit - > 5,000mg/kg

Acetic Acid Acute toxicity LD50 Oral – Rat – 3,310 mg/kg

LC50 Inhalation – Mouse – 1h - 5620ppm Remarks: Sense Organs and Special Senses (Nose, Eyes, Ears and Taste): Eyes: Conjunctive irritation. Eyes: Other. Blood: Other changes.

LD50 Dermal – Rabbit – 1,112 mg/kg

#### DMSO

**Skin corrosion/irritation** Skin – Rabbit – No skin irritation – 4h

Acetic Acid Skin corrosion/irritation Skin – Rabbit – Mild skin irritation – 24h

DMSO Serious eye damage/irritation Eyes – Rabbit – Mild eye irritation

Acetic Acid Serious eye damage/irritation Eyes – Rabbit – Corrosive to eyes.

**Respiratory or skin sensitisation** May cause sensitization by skin contact.

#### Germ cell mutagenicity

Genotoxicity in vitro – Mouse – lymphocyte Cytogenetic analysis Genotoxicity in vitro – Mouse – lymphocyte Mutation in mammalian somatic cells

Genotoxicity in vivo – Rat – Intraperitoneal Cytogenetic analysis

Genotoxicity in vivo - Mouse - Intraperitoneal





DNA damage

#### Carcinogenicity

Carcinogenicity – Rat – Oral Tumorigenic: Equivocal tumorigenic agent by RTECS criteria. Skin and Appendages: Others: Tumors.

Carcinogenicity – Mouse – Oral Tumorigic: Equivocal tumorigenic agent by RTECS criteria. Lukaemia skin and appendages: Other: Tumors.

IARC: No component of this product presents at levels greater than or equal to 0.1% is identified as probable, possible or confirmed human carcinogen by IARC.

#### **Reproductive toxicity**

Reproductive toxicity – Rat – Intraperitoneal Effects on fertility: Abortion

Reproductive toxicity – Rat – Intraperitoneal Effects on fertility: Post – implantation mortality (e.g. dead and/or resorbed implants per total number of implants).

Reproductive toxicity – Rat – Subcutaneous Effects on fertility: Post – implantation mortality (e.g. dead and/or resorbed implants per total number of implants). Effects on fertility: Litter size (e.g. # fetuses per litter; measured before birth).

#### Reproductive toxicity - Mouse - Oral

Effects on fertility: Pre-implantation mortality (e.g. reduction in number of implants per female; total number of implants per corpora lutea). Effects on Embryo or fetus: Fetoxicity (except death, e.g. stunted fetus). Specific developmental abnormalities: Musculoskeletal system.

#### Reproductive toxicity - Mouse - Intraperitoneal

Effects on embryo or fetus: Fetoxicity (except death, e.g. stunted fetus). Specific developmental abnormalities: Musculoskeletal system.

**STOT-single exposure** No data available

STOT-repeated exposure

No data available

Aspiration hazard. No data available

Potential Health Hazards	
Inhalation	Harmful if inhaled. Causes serious respiratory tract irritation.
Ingestion	Harmful if swallowed. Causes burns.
Skin	May be harmful if absorbed through skin. Causes skin burns.
Eyes	Causes eye irritation/ burns.
Aggravated Medical	
Condition	Avoid contact with DMSO solutions containing toxic materials or materials with unknown toxicological properties. Dimethyl sulfoxide is readily absorbed through the skin and may carry such materials into the body.

#### Signs and symptoms of exposure

Nausea, Fatigue, Headache. To the best of our knowledge, the chemical, physical and toxicological properties have not been thoroughly investigated.



Additional Information RTECS: PV6210000 RTECS: AF1225000

#### **SECTION 12. ECOLOGICAL INFORMATION**

12.1 Toxicity	
Toxicity to Fish 96h	LC50-Pimephales promelas (fathead minnow) – 34,000mg/l -
96h	LC50-Oncorhynchus mykiss (rainbow trout) – 34,000mg/l-
Toxicity to daphnia and other Aquatic invertebrates Toxicity to algae	EC50-Daphnia pulex (water fleas) – 27,500mg/l EC50-Lepomis macrochirus (bluegill) - >400,000mg/l-96h
Acetic Acid	
Toxicity to Fish	LC50 – Leuciscus idus (Golden Orfe) – 410.00mg/l – 48h LC50 – Cyprinus carpio (Carp) – 49.00mg/l – 48h LC50 – Pimephales promelas (Fathead minnow) – 79.00 -
88.00	Omg/I –
	960 LC50 – Lepomis macrochirus – 75mg/l – 96h
Toxicity to Daphnia and other aquatic invertebrates.	EC50 – Daphnia magna (Water flea) – 65.00mg/l – 48h
12.2 Porsistones and degradabil	ity

#### **12.2 Persistence and degradability** Biodegradability

Remarks: Expected to be biodegradable.

#### 12.3 Bioaccumulative potential

No data available

#### 12.4. Mobility in soil

No data available

## 12.5. Results of PBT and vPvB assessment

No data available

## 12.6. Other adverse effects

Biochemical Oxygen Demand (BOD) - 880mg/g

# **SECTION 13. DISPOSAL CONSIDERATIONS**

#### 13.1 Waste treatment methods

This combustible material may be burned in a chemical incinerator equipped with an afterburner and scrubber or to be disposed of by a licensed professional waste disposal company.

#### Contaminated packaging

Dispose of as the unused product.

# **SECTION 14. TRANSPORT INFORMATION**

#### 14.1 UN Number

DMSO Acetic Acid	ADR/RID: - ADR/RID: 278	9	IMDG: - IMDG: 2789	IATA: - IATA: 2789
14.2 UN Prop	er Shipping N	ame	) and a	
DIVISO	IMDG:	Not Dangerous G	Goods	
Acetic Acid	IATA: ADR/RID: IMDG: IATA:	Not Dangerous G ACETIC ACID, G ACETIC ACID, G Acetic Acid, glaci	Goods LACIAL LACIAL al	
14.3 Transpo	rt hazard clas	s (es)		
DMSO Acetic Acid	ADR/RID: - ADR/RID: 8 (3	3)	IMDG: - IMDG: 8 (3)	IATA: - IATA: 8 (3)
14.4 Packing	group			
DMSO Acetic Acid	ADR/RID: - ADR/RID: II		IMDG: - IMDG: II	IATA: - IATA: II
14.5 Environr	<b>nental hazard</b> ADR/RID: No	S	IMDG Marine pollutant: No	IATA: No
14.6 Special	14.6 Special precautions for user			

No data available

# **SECTION 15. REGULATORY INFORMATION**

This safety data sheet complies with the requirements of Regulation (EC) No. 1907/2006

# 15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture

No data available

## **15.2 Chemical Safety Assessment**

No data available

Please note that the label elements that used to go in Section 15 are now in Section 2.

# **SECTION 16. OTHER INFORMATION**

The advice offered is derived from the current available information on the hazardous materials in this product and it component(s). Consideration has been made regarding the quantities offered in the pre dispensed container. The advice offered is, therefore not all inclusive nor should it be taken as the descriptive of the compound generally.

# SAFETY DATA SHEET

Version: 2.0 Date written: 12<sup>th</sup> January 2015 Date reviewed: 27<sup>th</sup> April 2018

# 1. IDENTIFICATION OF THE SUBSTANCE/PREPARATION AND OF THE COMPANY / UNDERTAKING

Product Name	Procainamide Dye
Product Catalogue Name	LT-PROC-01, LT-PROC-96
CAS-No.	614-39-1
Company:	Ludger Ltd Culham Science Centre Abingdon Oxfordshire OX14 3EB
Telephone: Emergency Telephone: Email:	01865 408554 01865 408554 info@ludger.com

# 2. HAZARDS IDENTIFICATION

#### 2.1 Classification of the substance or mixture Classification according to Regulation (EC) No 1272/2008 [EU-GHS/CLP] Acute toxicity, Oral (Category 4) Skin irritation (Category 2) Eye irritation (Category 2) Specific target organ toxicity - single exposure (Category 3)

## 2.2 Label elements



Signal Word: Warning

## Hazard Statement(s)

H302 Harmful if swallowed.H315 Causes skin irritation.H319 Causes serious eye irritation.H335 May cause respiratory irritation.

## **Precautionary Statement(s)**

P261 Avoid breathing dust/ fume/gas/mist/vapours/spray P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do so. Continue rinsing.

#### 2.3 Other hazard information:



No supplemental hazard statements.

#### **3. COMPOSITION/INFORMATION ON INGREDIENTS**

### 3.1 Substances

Synonyms:	Procainamide hydrochloride; 4-Amino-N-(2-diethylaminoethyl)benzamide
	hydrochloride; 4-Aminobenzoic acid 2-diethylaminoethylamide
Formula:	
Molecular Weight:	Procainamide hydrochloride: 271.79 g/mol

Component		Concentration
Name	Procainamide Dye	-
CAS-No.	614-39-1	
EC-No.	210-381-7	

## **SECTION 4. FIRST AID MEASURES**

#### 4.1 Description of first aid measures

#### **General Advice**

Consult a physician if exposure causes ill effects and if in any doubt. Show this safety data sheet to the physician/ first responder in attendance.

#### If Ingested

Never give anything by mouth to an unconscious person. Rinse mouth with water. Consult a physician.

#### If skin is exposed

Wash off with soap and plenty of water. Consult a physician.

#### If eyes are exposed

Rinse thoroughly with plenty of water for at least 15 minutes and consult a physician.

#### If inhaled

If breathed in, move person into fresh air. If not breathing, give artificial respiration. Consult a physician.

#### 4.2 Most important symptoms and effects, both acute and delayed

The most important known symptoms and effects are described in the labelling (see section 2.2) and/or in

# section 11

# 4.3 Indication of immediate medical attention and special treatment needed

No data available.

## 5. FIRE-FIGHTING MEASURES

#### 5.1 Extinguishing media

Use water spray, alcohol-resistant foam, dry chemical or carbon dioxide.

#### 5.2 Special hazards arising from the substance or mixture

Carbon oxides, nitrogen oxides (NOx), Hydrogen chloride gas.

#### 5.3 Advice for firefighters

Wear self contained breathing apparatus for fire fighting if necessary.

## 6. ACCIDENTAL RELEASE MEASURES

#### 6.1 Personal precautions, protective equipment and emergency procedures

Use personal protective equipment. Avoid dust formation. Avoid breathing vapours, mist or gas. Ensure

adequate ventilation. Evacuate personnel to safe areas. Avoid breathing dust.

#### **6.4 Environmental Precautions**

Do not let product enter drains.

#### 6.5 Methods and material for containment and cleaning up

Pick up and arrange disposal without creating dust. Sweep up and shovel. Keep in suitable, closed containers for disposal.

#### 6.4 Reference to other sections

Please refer to section 13.

# 7. HANDLING AND STORAGE

#### 7.1 Precautions for safe handling

Avoid contact with skin and eyes. Avoid formation of dust and aerosols. Provide appropriate exhaust ventilation at places where dust is formed.Normal measures for

preventive fire

# protection.

#### 7.2 Conditions for safe storage, including any incompatibilities

Store in a cool, dark place. Keep the container tightly closed in a dry well ventilated place.

#### 7.3 Specific end uses

A part from the uses mentioned in section 1.2 no other specific uses are stipulated.

#### 8. EXPOSURE CONTROLS/PERSONAL PROTECTION

#### 8.1 Control parameters

Contains no substances with occupational exposure limit values.

#### 8.3 Exposure controls

#### Appropriate engineering controls

Handle in accordance with good industrial hygiene and safety practice. Wash hands before breaks and at the end of workday.

#### **Personal Protective Equipment**

#### Eye / face protection

Safety glasses with side-shields conforming to EN166 Use equipment for eye protection tested and approved under appropriate government standards such as NIOSH (US) or EN 166(EU).

#### Skin protection



Handle with gloves. Gloves must be inspected prior to use. Use proper glove removal technique (without touching glove's outer surface) to avoid skin contact with this product. Dispose of contaminated gloves after use in accordance with applicable laws and good laboratory practices. Wash and dry hands.

The selected protective gloves have to satisfy the specifications of EU Directive 89/686/EEC and the standard EN 374 derived from it.

#### **Body Protection**

Complete suit protecting against chemicals. The type of protective equipment must be selected according to the concentration and amount of the dangerous substance at the specific workplace.

#### **Respiratory protection**

For nuisance exposures use type P95 (US) or type P1 (EU EN 143) particle respirator.For higher level protection use type OV/AG/P99 (US) or type ABEK-P2 (EU EN 143) respirator cartridges. Use respirators and components tested and approved under appropriate government standards such as NIOSH (US) or CEN (EU).

# 9. PHYSICAL AND CHEMICAL PROPERTIES

#### 9.1 Information on basic physical and chemical properties

Appearance Odour Odour threshold bН Freezing/Melting Point Initial boiling point and boiling range Flash Point **Evaporation** rate Flammability Upper/lower flammability or explosive limits Vapour Pressure, Pa at temperature degree C **Relative Densitv** Solubility in water and solvents (mg/l) Partition coefficient Autoignition temperature Decomposition temperature Viscosity Explosive properties Oxidising properties

Form: solid No data available No data available No data available Melting point/range: 167 - 169 °C - lit. No data available No data available

#### 9.2 Other information

No data available

## **10. STABILITY AND REACTIVITY**

10.1 Reactivity

No data available

#### **10.2 Chemical stability**

Stable under recommended storage conditions. **10.3 Possibility of hazardous reactions** No data available



#### **10.4 Conditions to Avoid**

No data available

#### **10.5 Incompatible materials**

Strong oxidizing agents

# 10.6 Hazardous decomposition products

Other decomposition products - No data available

# **11. TOXICOLOGICAL INFORMATION**

#### 11.1 Information on toxicological effects

Acute toxicity R LD50 Oral - rat - 1,509 mg/kg

#### Skin corrosion/irritation No data available

Serious eye damage/irritation No data available

Respiratory or skin sensitisation No data available

Germ cell mutagenicity No data available

#### Carcinogenicity

No component of this product present at levels greater than or equal to 0.1% is identified as probable, possible or confirmed human carcinogen by IARC.

**Reproductive toxicity** No data available

# STOT-single exposure

Inhalation - May cause respiratory irritation. **STOT-repeated exposure** No data available

#### Aspiration hazard.

No data available

#### **Potential Health Hazards**

Inhalation	May be harmful if inhaled. Causes respiratory tract irritation.
Ingestion	Harmful if swallowed.
Skin	May be harmful if absorbed through the skin. Causes skin irritation.
Eyes	Causes serious eye irritation.

## Signs and symptoms of exposure

To the best of our knowledge the chemical, physical and toxicological properties have not been thoroughly investigated.



# **12. ECOLOGICAL INFORMATION**

#### 12.1 Toxicity

No data available

### **12.2 Persistence and degradability**

No data available

### 12.3 Bioaccumulative potential

No data available

#### 12.4. Mobility in soil

No data available

#### 12.5. Results of PBT and vPvB assessment

PBT/vPvB assessment not available as chemical safety assessment not required/not conducted

#### 12.6. Other adverse effects

No data available

# **13. DISPOSAL CONSIDERATIONS**

#### 13.1 Waste treatment methods

Contact a licensed waste disposal service to collect/dispose of any waste material.

#### Contaminated packaging

Treat as an unopened/ unused product.

#### **14. TRANSPORT INFORMATION** 14.1 UN Number ADR/RID: -IMDG: -IATA: -14.2 UN Proper Shipping Name Not dangerous goods ADR/RID: IMDG: Not dangerous goods IATA: Not dangerous goods 14.3 Transport hazard class(es) ADR/RID: -IMDG: -IATA: -14.4 Packing group IATA: -ADR/RID: -IMDG: -14.5 Environmental hazards no IMDG Marine pollutant: no IATA: no ADR/RID: 14.6 Special precautions for user No data available

# **15. REGULATORY INFORMATION**

This safety datasheet complies with the requirements of Regulation (EC) No. 1907/2006.



# 15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture

No data available

#### **15.2 Chemical Safety Assessment**

For this product a chemical safety assessment was not carried out.

Note that the label elements, the Risk and Safety phrases (now Hazard and Precautionary statements) that used to be in Section 15 are now in Section 2.

## **16. OTHER INFORMATION**

The advice offered is derived from the current available information on the hazardous materials in this product and it component(s). Consideration has been made regarding the quantities offered in the pre dispensed container. The advice offered is, therefore not all inclusive nor should it be taken as the descriptive of the compound generally.

# SAFETY DATA SHEET

Version: 3.0 Date written: 1<sup>st</sup> February 2012 Date reviewed: 27<sup>th</sup> April 2018

# SECTION 1. IDENTIFICATION OF THE SUBSTANCE/PREPARATION AND OF THE COMPANY/ UNDERTAKING

Product Name	sodium cyanoborohydride		
Product Catalogue Name	LT-CYANOB-03, LT-CYANOB-05, LT-CYANOB-96		
CAS-No.	25895-60-7		
Company:	Ludger Ltd Culham Science Centre Abingdon Oxford OX14 3EB		
Telephone: Emergency Telephone: Email:	01865 408554 01865 408554 info@ludger.com		

## **SECTION 2. HAZARDS IDENTIFICATION**

#### 2.1 Classification of the substance or mixture Classification according to the Regulation (EC) No. 1272/2008 [EU-GHS/CLP]

Flammable solids (category 1) Substances, which in contact with water emit flammable gases (category 1) Acute toxicity, Oral (category 2) Acute toxicity, Inhalation (category 2) Acute toxicity, Dermal (category 2) Skin corrosion (category 1B) Acute aquatic toxicity (category 1) Chronic aquatic toxicity (category 1)

## 2.2 Label elements



Signal Word: DANGER

Hazard Statement(s)	Flammable solid
11220	Tranninable Solic
H260	In contact with water releases flammable gases which may ignite
spontaneously.	
H300	Fatal if swallowed.
H310	Fatal in contact with skin.
H314	Causes severe burns and eye damage.
H330	Fatal if inhaled.
H410	Very toxic to aquatic life with long lasting effects.

## Precautionary Statement(s)

P210 Keep away from heat/sparks/open flames/hot surfaces. No smoking.



P223 and	Keep away from any possible contact with water, because of violent reaction
	possible flash fire.
P260	Do not breathe dust/ fume/ gas/ mist/ vapours/spray.
P264	Wash hands thoroughly after handling.
P273	Avoid release to the environment.
P280	Wear protective gloves/ protective clothing/ eye protection/ face protection.
P284	Wear respiratory protection.

#### 2.3 Other hazard information (EU) EUH032 Contact v

Contact with acids liberates very toxic gas.

## **SECTION 3. COMPOSITION/INFORMATION ON INGREDIENTS**

#### 3.1 Substances

Synonyms:	Sodium Cyanotrihydridoborate
Formula:	CH3BNNa
Molecular weight	62.84 g/mol

Component		Concentration
Sodium Cy	anotrihydroborate	
CAS-No.	25895-60-7	-
EC-No.	247-317-2	

# **SECTION 4. FIRST AID MEASURES**

#### 4.1 Description of first aid measures

#### **General Advice**

Consult a physician if exposure causes ill effects and if in any doubt. Show this safety data sheet to the doctor/ first responder in attendance.

#### If Ingested

Do NOT induce vomiting. Rinse mouth well with water, unless person(s) is unconscious.

#### If skin is exposed

Remove contaminated clothing/shoes immediately. Wash affected area(s) with water and soap.

#### If eyes are exposed

Wash eye(s) with plenty of water for at least 15 minutes, if unsure seek medical advice.

#### If inhaled

Move into a source of fresh air, if not breathing give artificial respiration.

#### 4.2 Most important symptoms and effects, both acute and delayed

Burning sensation, coughing, wheezing, laryngitis, shortness of breath, head ache, nausea, vomiting.

To the best of our knowledge the chemical, physical and toxicological properties have not been thoroughly investigated. Absorption into the body leads to the formation of methemoglobin which in sufficient concentration causes cyanosis.

Onset may be delayed two to four hours or longer.

## 4.3 Indication of immediate medical attention and special treatment needed

No Data available

# **SECTION 5. FIRE-FIGHTING MEASURES**



#### 5.1 Extinguishing media

Dry powder

#### 5.2 Special hazards arising from the substance or mixture

Carbon oxides, Nitrogen oxides, Hydrogen cyanide (Hydrocyanic acid), Borane/boron oxides.

#### 5.3 Advice for firefighters

Wear self contained breathing apparatus for fire fighting if necessary.

### **SECTION 6. ACCIDENTAL RELEASE MEASURES**

#### 6.1 Personal precautions, protective equipment and emergency procedures.

Wear respiratory protection; gently sweep up to avoid dust formation. Avoid breathing vapours, mist or gas. Ensure adequate ventilation, remove all sources of ignition. Evacuate personnel to a safe area; avoid breathing in dust/gas or mist.

For personal protection see section 8.

#### **6.2 Environmental Precautions**

Prevent further leakage or spillage if safe to do so. Do not let the chemical enter the drainage system and further discharge into the environment must be avoided.

#### 6.3 Methods and material for containment and cleaning up

Contain the spill with matting if necessary and then collect using either an electrically protected vacuum cleaner or by damp brushing (not wet) and putting the collected waste into a secure dry container, do not flush with water. Dispose according to local regulations.

#### 6.4 Reference to other sections

For disposal regulations see section 13.

#### **SECTION 7. HANDLING AND STORAGE**

#### 7.1 Precautions for safe handling

Avoid contact with skin and eyes, avoid formation of dust when handling. Provide appropriate exhaust ventilation in work areas where dust could be formed. Keep away from sources of ignition (No Smoking) and take measures to prevent the build up of electrostatic charge.

#### 7.2 Conditions for safe storage, including any incompatibilities

Store in a cool place. Keep the container tightly closed in a dry and well ventilated place. Never allow the product to get into contact with water during storage as it is moisture sensitive. Do not store near acids. Handle and open the container with care. Hygroscopic. Handle when open under an inert gas.

#### 7.3 Specific end uses

No data available

# SECTION 8. EXPOSURE CONTROLS/PERSONAL PROTECTION

#### 8.1 Control parameters

Contains no substances with occupational exposure limit values.

Component	CAS No.	Value Form of exposure	Control parameters	Basis
sodium	25895-60-7	TWA	5 mg/m3	UK. EH40 WEL -



cyanoborohydride				Workplace
				Exposure Limits
	Remarks	Can be absorbed through skin. The assigned substances are those for which there are concerns that dermal absorption will lead to systemic toxicity. Where no specific short-term exposure limit is listed, a figure threetimes the long-term exposure should be used		

#### 8.2 Exposure controls Appropriate engineering controls

Avoid contact with skin, eyes and clothing. Wash hands before breaks and immediately after handling the product.

# Personal protective equipment

## Eye/Face protection

Face shield and safety glasses to be worn following good laboratory practice. Eye protection should be tested and approved under appropriate government standards such as EN 166 (EU) or NIOSH (US).

# Skin protection

Handle with gloves at all times following good laboratory practice. Gloves must be inspected before use and to be removed in proper glove removal technique (without touching the gloves outer surface) to avoid skin contact. Dispose of contaminated gloves after use as contaminated waste, in accordance with local regulations. Wash and dry hands.

Gloves to be within the specifications of EU directive 89/686/EEC and the standard EN 374 derived from it.

Full contact Material: Nitrile rubber Minimum layer thickness: 0.11 mm Break through time: 480 min Material tested:Dermatril® (KCL 740 / Aldrich Z677272, Size M)

Splash contact Material: Nitrile rubber Minimum layer thickness: 0.11 mm Break through time: 480 min Material tested:Dermatril® (KCL 740 / Aldrich Z677272, Size M) data source: KCL GmbH, D-36124 Eichenzell, phone +49 (0)6659 87300, e-mail sales@kcl.de, test method: EN374

If used in solution, or mixed with other substances, and under conditions which differ from EN 374, contact the supplier of the CE approved gloves. This recommendation is advisory only and must be evaluated by an industrial hygienist and safety officer familiar with the specific situation of anticipated use by our customers. It should not be construed as offering an approval for any specific use scenario.

## **Body protection**

Complete suit protecting against chemicals, Flame retardant antistatic protective clothing. The type of protective equipment must be selected according to the concentration and amount of the dangerous substance at the specific workplace.

## **Respiratory protection**

Where risk assessment shows air-purifying respirators are appropriate use a full-face particle



respirator type N100 (US) or type P3 (EN 143) respirator cartridges as a backup to engineering controls. If the respirator is the sole means of protection, use a full-face supplied air respirator. Use respirators and components tested and approved under appropriate government standards such as NIOSH (US) or CEN (EU).

#### **Control of environmental exposure**

Prevent further leakage or spillage if safe to do so. Do not let product enter drains. Discharge into the environment must be avoided.

# SECTION 9. PHYSICAL AND CHEMICAL PROPERTIES

#### 9.1 Information on basic physical and chemical properties

Appearance	Form: Powder
	Colour: Beige
Odour	No data available
Odour threshold	No data available
рН	No data available
Freezing/Melting Point	Melting point/range :> 242° C
Initial boiling point and boiling range	No data available
Flash Point	No data available
Evaporation rate	No data available
Flammability (solids and gases)	The substance or mixture is a flammable solid
	with the subcategory 1.
Upper/lower flammability or explosive limits	No data available
Vapour Pressure	No data available
Relative Density	No data available
Solubility in water and solvents	No data available
Partition coefficient	No data available
Autoignition temperature	No data available
Decomposition temperature	No data available
Viscosity	No data available
Explosive properties	No data available
Oxidising properties	No data available

#### 9.2 Other information

Oxidising properties

None available

## SECTION 10. STABILITY AND REACTIVITY

#### **10.1 Reactivity**

No data available

#### **10.2 Chemical stability**

No data available

#### **10.3 Possibility of hazardous reactions**

Reacts violently with water.

#### **10.4 Conditions to Avoid**

Do not allow water to enter container because of violent reaction. Heat, flames and sparks. Extremes of temperature and direct sunlight. Exposure to sunlight.

#### **10.5 Incompatible materials**

Do not store near acids, oxidising agents.



#### **10.6 Hazardous decomposition products**

Hazardous decomposition products formed under fire conditions. - Carbon oxides, Nitrogen oxides (NOx), Hydrogen cyanide (hydrocyanic acid), Borane/boron oxides Reacts with water to form: - Hydrogen gas Other decomposition products - No data available Hazardous decomposition products formed under fire conditions. - Carbon oxides, Borane/boron oxides, Sodium oxides

# SECTION 11. TOXICOLOGICAL INFORMATION

#### 11.1 Information on toxicological effects Acute toxicity

No Data available

#### Skin corrosion/irritation

No data available

Serious eye damage/irritation

No data available

**Respiratory or skin sensitisation** No data available

Germ cell mutagenicity No data available

#### Carcinogenicity

**IARC:** No component of this product present at levels greater than or equal to 0.1% is identified as probable, possible or conformed human carcinogen by IARC.

#### **Reproductive toxicity**

No data available

#### STOT (specific target organ toxicity) -single exposure

No data available

# STOT (specific target organ toxicity) -repeated exposure

No data available

#### Aspiration hazard.

No data available

#### **Potential health effects**

Inhalation	May be fatal if inhaled. Material is extremely destructive to the tissue
of the	
	mucus membranes and upper respiratory tract.
Ingestion	May be fatal if swallowed. Causes burns.
Skin	May be fatal if absorbed through the skin. Causes skin burns.
Eyes	Causes burns to the eyes.

#### Signs and symptoms of Exposure

Burning sensation, coughing, wheezing, laryngitis, shortness of breath, head ache, nausea, vomiting. To the best of our knowledge the chemical, physical and toxicological properties have not been thoroughly investigated. Absorption into the body leads to the formation of methemoglobin which in sufficient concentration causes cyanosis. Onset may be delayed two to four hours or longer.



# **SECTION 12. ECOLOGICAL INFORMATION**

**12.1 Toxicity** No data available

**12.2 Persistence and degradability** No data available

**12.3 Bioaccumulative potential** No data available

12.4. Mobility in soil

No data available

**12.5. Results of PBT and vPvB assessment** No data available

#### 12.6. Other adverse effects

Very toxic to aquatic life with long lasting effects.

#### **SECTION 13. DISPOSAL CONSIDERATIONS**

#### 13.1 Waste treatment methods

Product or/and collected waste from spillage.

Burn in a chemical incinerator equipped with an afterburner and scrubber but take extra precautions when igniting as this material is highly flammable. Or to contact a licensed disposal company and arrange disposal, inform the company of the nature of the waste.

#### Contaminated packaging

Dispose of as the unused product, with a licensed disposal company.

## **SECTION 14. TRANSPORT INFORMATION**

#### 14.1 UN Number ADR/RID: 3179 IMDG: 3179 IATA: 3179 14.2 UN Proper Shipping Name FLAMMABLE SOILD, TOXIC, INORGANIC, N.O.S. (Sodium cyanotrihydroborate) ADR/RID: IMDG: FLAMMABLE SOILD, TOXIC, INORGANIC, N.O.S. (Sodium cyanotrihydroborate) IATA: Flammable solid, toxic, inorganic, n.o.s. (Sodium cyanotrihydroborate) 14.3 Transport hazard class (es) ADR/RID: 4.1 (6.1) IMDG: 4.1 (6.1) IATA: 4.1 (6.1) 14.4 Packing group ADR/RID: II IMDG: II IATA: II 14.5 Environmental hazards ADR/RID: No IMDG Marine pollutant: No IATA: No 14.6 Special precautions for user

No data available

# **SECTION 15. REGULATORY INFORMATION**



# 15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture

No data available

# 15.2 Chemical Safety Assessment

No data available

# **SECTION 16. OTHER INFORMATION**

The advice offered is derived from the current available information on the hazardous materials in this product and it component(s). Consideration has been made regarding the quantities offered in the pre dispensed container. The advice offered is, therefore not all inclusive nor should it be taken as the descriptive of the compound generally.