



Permethylation of Glycans

LudgerTag™ technology to enable
rapid, reliable, high-throughput (HT)
MALDI-TOF-MS analysis

 Ludger



Glycan Permethylation



Permethylation involves the addition of methyl groups (CH_3) to all of the hydroxyl and *N*-acetyl groups, and also methyl esterifies the carboxy function on the sialic acid.


Why Permethylate?

(released N- and O-glycans)


Improves and enhances Ionization efficiency of
glycans on mass spec (MALDI-MS, ESI-MS)

When compared to non-derivatized oligosaccharides

Increased glycan
hydrophobicity enables
LC-MS analysis



↓ CH₃I, NaOH




Easier determination of
branching and glycosidic
linkage positions

Relative and absolute
quantitation can be performed
*By introducing isotope
labelled internal standards*

Enables detection of both
neutral and acidic glycans
in positive ion mode using
MALDI-TOF-MS

Stabilizes the labile
sialic acid moieties

Conventional Permethylation Methods are Labour Intensive



Picture taken from: Akihiko Kameyama, (2014). GlycoPOD <http://jcgddb.jp/GlycoPOD>. Web.15,8,2014

Slurry method

A slurry of sodium hydroxide in dimethyl sulfoxide needs to be crushed and prepared as shown in the illustration and this slurry needs to be constantly vortexed before adding it to individual sample vials.

This approach is **slow, low throughput and labour intensive**.

Solid phase methods

Solid phase permethylation techniques have been developed more recently however, many of these **techniques are still labour intensive and repetitive** which are **not practical** for large sample numbers (e.g. for the characterization of biopharmaceuticals).

Reference:

- (1) Pilsoo Kang, Y. M. and M. V. N. *RAPID Commun. MASS Spectrom.* 2008, 22, 721–734.
- (2) Jeong, H.-J.; Kim, Y.-G.; Yang, Y.-H.; Kim, B.-G. *Anal. Chem.* 2012, 84 (7), 3453–3460.
- (3) Gao, X.; Zhang, L.; Zhang, W.; Zhao, L. *Analyst* 2015, 140 (5), 1566–1571.

RAPID COMMUNICATIONS IN MASS SPECTROMETRY
Rapid Commun. Mass Spectrom. 2008; 22: 721–734
Published online in Wiley InterScience (www.interscience.wiley.com) DOI: 10.1002/rcm.3395

RCM

High-throughput solid-phase permethylation of glycans prior to mass spectrometry

Pilsoo Kang, Yehia Mechref* and Milos V. Novotny*

Department of Chemistry, Indiana University, Bloomington, IN 47405, USA

Received 6 November 2007; Accepted 17 December 2007

Permethylation of glycans prior to their mass spectrometric determination has now become a time-honored methodology in glycoconjugate analysis due to the advantage of a simultaneous analysis of neutral and acidic glycans as well as enhanced sensitivity and easier tandem mass spectrometry interpretation. While the different solvent extraction-based versions of this method often suffice in different structural studies, they are generally less satisfactory in the quantitative determinations aiming at minor quantities of the analyzed materials. To overcome these difficulties, we recently introduced a solid-phase capillary permethylation technique (Kang *et al.*, *Rapid Commun. Mass Spectrom.* 2005; 19: 3421) for microscale determination. Here, we describe a very useful high-throughput extension of the solid-phase methodology utilizing spin columns packed with sodium hydroxide beads. This procedure has been thoroughly optimized to match the analytical

Ludger Glycan Permethylation Kit (LT-PERMET-96)

Simple to use

Kit format is less labour intensive than conventional in-solution permethylation

Data is comparable to gold standard HILIC UHPLC data
(2AB / Procainamide labelling)

A microplate based 96-well plate
Format of kit is convenient as it is scalable between 1 to 96 samples.



Reliable and validated according to EMA ICH Q2 (R1) guidelines

Compatible with MALDI-TOF-MS and LC-MS
Absolute quantitation can be achieved by introducing Isotope labelled internal standards.

Method can be automated for HT, rapid sample prep
Workflow can be adapted to a liquid handling robot

Suitable for N- and O-glycan analysis and also aids linkage analysis
Easier determination of branching and glycosidic linkage analysis

Components of the LT-PERMET-96 Kit

1. Derivatisation



Addition of 300µL DMSO to LT-PERMET-96 plate (15 min incubation)

Dimethyl sulfoxide (DMSO)
LT-PERMET-DMSO-96



Addition of 55µL methyl iodide (1 hour incubation)

Methyl Iodide (MeI)
LT-PERMET-MeI-96

2. Liquid Liquid Extraction



Addition of 400µL DCM and 1ml water

Dichloromethane (DCM)
LT-PERMET-DCM-96




Performing LLE to render the pH neutral prior to MALDI-TOF-MS measurement

A black and white photograph showing a hand holding a microcentrifuge tube. The tube is held between the thumb and index finger, with the middle finger supporting it from the side. The tube is partially filled with a clear liquid. The background is dark and out of focus.


Permethylatation Workflow
using the LT-PERMET-96 kit

Workflow - LudgerTag Permethylation for N-glycan Profiling and Identification




Native Glycans	Fucose	Galactose	Mannose	N-Acetylglucosamine	N-Acetylneuraminic acid
Permethylated Glycans	Fucose	Galactose	Mannose	N-Acetylglucosamine	N-Acetylneuraminic acid

Workflow - LudgerTag Permethylation for O-glycan Profiling and Identification



Native Glycans



Permethylated Glycans



Recommended Components for Workflow – LudgerTag Permethylation for N- and O-glycan Profiling

Product	Product code	N- glycans	O- glycans
PNGase F N-glycan release	E-PNG-01	●	
Hydrazinolysis or Orela O-glycan release kit	LL-HYDRAZ-A2* Technical guide		●
	LL-ORELA-A2 Technical guide		●
Enrichment plate	LC-PERMET-96 Technical guide	●	
Cation exchange clean-up cartridges	LC-CEX-A6 Technical guide		●
Permethylation kit	LT-PERMET-96 Technical guide	●	●


*Note: although hydrazinolysis can also be used for N glycan release, we recommend PNGase F release for the workflow




LT-PERMET-96:
Manual Workflow

LT-PERMET-96: Manual Procedure


Addition of DMSO, Mel and Incubation




Addition of DCM and water




Liquid-liquid extraction



Ready for Analysis





LT-PERMET-96: Automated Workflow

The system can be adapted to a liquid handling robot,
enabling reliable High Throughput (HT) Studies

LT-PERMET-96: Automated HT Workflow

Typical timeline for 96 samples




2h
(16h incubation 37C°)

2.5h
(Plus drying time)

1h
(Plus 1.5h incubation)

3h
(Plus drying time)

~1min/sample = 1.6h
for data acquisition








1. Automated
PNGase-F release

2. Automated
HILIC-SPE enrichment

3a. Automated
Permethylat

3b. Automated Liquid-
Liquid Extraction (LLE)

4. Data acquisition
(MALDI-TOF-MS)





Examples of how
we're using the
LT-PERMET-96 system
at Ludger

Permethylated and 2-AB Labelled Human IgG N-glycan Profiles are Comparable




Figure 1. Comparison of PNGase F released and purified human IgG N-glycans analyzed with orthogonal methods (A) MALDI-TOF-MS spectrum of permethylated human-IgG N-glycans (B) 2-AB labelled HILIC UHPLC chromatogram.

▲ Fucose ● Galactose ● Mannose ■ *N*-Acetylglucosamine ◆ *N*-Acetylneuraminic acid

- We analysed human IgG N-glycans using the automated HT permethylation method and compared the data to those obtained from UHPLC analysis by fluorescence detection as shown in Figure 1 Section A and B.
- Glycan signals were integrated, normalized and the relative intensities and standard deviation was calculated for the 13 major N-glycan peaks.
- The analysis confirmed that peaks with higher relative intensities (above 4%) showed good correlation between the two methods.
- Therefore we conclude that the HT permethylation technique is comparable to UHPLC results and that it gives a reliable overview of the glycosylation profile in a short timespan.

Permethylation Stabilises Fragile Sialic Acids




Figure 2. Comparison of the glycosylation profiles of the A3G3S3 N-glycan standard analyzed after sample preparation using the liquid handling robot. (A) MALDI-TOF-MS spectrum after permethylation, (B) HILIC UHPLC chromatogram with fluorescence detection after procainamide labelling. (C) Histogram comparing the relative peak intensities of triantennary, disialylated structures (H6N5S2) and triantennary, trisialylated structures (H6N5S3) after triplicate analysis. The histogram shows comparable relative signal intensities between MALDI-TOF-MS and HILIC UHPLC analysis. The error bars depict standard deviation.


- The relative intensities of triantennary, disialylated structures (H6N5S2) and triantennary, trisialylated structures (H6N5S3) from the A3G3S3 glycan standard were determined by MALDI-TOF-MS after automated HT permethylation.
- This data was then compared to the ratios obtained after procainamide labeling followed by HILIC UHPLC with fluorescence detection as shown in Figure 2.
- Triplicate analysis and relative quantitation was performed for both fluorescent labeling and permethylation.
- The analysis confirmed that the MALDI-TOF-MS data from the sialylated N-glycan standard gave similar and comparable results to that of the UHPLC data.

How to Start Using the Ludger Glycan Permethylation Technology



1. Submit your samples for automated HT permethylation

As part of our glycoprofiling services, we can perform sample preparation and analysis for you in our labs



2. Method transfer

We can transfer the glycan permethylation methods to your lab and provide technical support



3. Use our permethylation kit in your own lab

*Contact us for a quotation and place your order
Catalogue # LT-PERMET-96*

Contact Us

If you have
technical questions



CLICK
to contact
Archana

Archana Shubhakar
Senior Scientist
archana.shubhakar@ludger.com

To request a quotation:
For services, method transfer
or LT-PERMET-96 kit



CLICK
to contact
Sales

Sales Team
Quotations: info@ludger.com
Orders: sales@ludger.com