


Applications of Ludger GPEP IgG Standard


1) Process Control



2) System Suitability Standard

Glycan structure	Glycan Name	Peak ID
	A1F	16, 18, 19
	FA2	1, 3, 7
	FA2B	2, 10
	FA2G1	4, 6, 8, 9, 12, 14
	FA2G1B	5, 13
	NA2F	11, 15, 17

Liquid Chromatography-Mass Spectrometry (LC-MS)



Matrix-Assisted Laser Desorption/Ionization (MALDI)

