





LudgerTag Procainamide N-Glycan Analysis Workflow

 Glycoprotein	 Released N-Glycans + Proteins	 Labelled N-Glycans + Reagents + Reaction + Byproducts	 Enriched Labelled N-Glycans
Release (15 min)	Labelling (1 hour)	Clean Up (1 hour)	Analysis
PNGaseF release kit E-PNG01 E-rPNG01 LZ-rPNGaseF-kit	Procainamide glycan labelling kit LT-KPROC-24 LT-KPROC-96 LT-KPROC-VP24	LudgerClean procainamide plate LC-PROC-96	(U)HPLC, LC-MS or MALDI-MS
Standards & Controls (run with your samples)			
Human IgG Glycoprotein Standard GCP-IGG-100U			Procainamide-labelled glucose homopolymer ladder CPROC-GHP-30